

INVOS® Cerebral/Somatic Oximeter Quick Reference Guide for Pediatric Use




SOMANETICS® CEREBRAL/SOMATIC
INVOS OXIMETER
REFLECTING THE COLOR OF LIFE®


Key Code



Operating Basics



Science and Interventions



Case Examples



These guidelines are not designed to replace clinical judgment or individual patient needs. For complete instructions, indications for use, warnings and precautions, see the Operations Manual.

Property of Somanetics Corporation. Duplication of this guide is prohibited without written authorization from Somanetics Corporation.

Description of Symbols

OXIMETER ENCLOSURE



Blue LED – indicating AC Power Supply and charge ON



Caution; consult accompanying documents



Digital Output



Green LED – indicating power ON



HOME



ON/OFF



USB Port



VGA Port

DISPLAY



Alarm



Baseline



Battery



Signal Strength Indicator



Silence Alarm

OTHER LABELS



Consult operating instructions




Do not reuse

The noninvasive INVOS Cerebral / Somatic Oximeter is intended for use as an adjunct trend monitor of regional hemoglobin oxygen saturation of blood in the brain of an individual. It is also intended for use as an adjunct trend monitor of hemoglobin oxygen saturation of blood in a region of skeletal muscle tissue beneath the sensor in infants, children or adults at risk for reduced-flow or no-flow ischemic states. The prospective clinical value of data from the INVOS System has not been demonstrated in disease states. The INVOS System should not be used as the sole basis for diagnosis or therapy.



Initial Setup

1. Plug in INVOS Monitor to power source.

THE BLUE  LED ILLUMINATES TO INDICATE THE POWER IS CONNECTED AND THE BATTERY IS CHARGING. IF BLUE LIGHT IS NOT ON, TURN BACK PANEL MASTER (MAINS) POWER SWITCH ON.

2. Connect Preamplifier(s) to INVOS Monitor. Align red dot on the silver connector with red dot on side panel input connection. Insert cable connector into side panel connection while keeping the dots aligned. Be sure to fully insert the connector until it locks.



3. Connect the Reusable Sensor Cable Connectors to Preamplifier(s). Use color-coding.




Initial Setup (continued)

4. Attach Sensors to Reusable Sensor Cables. Hold blue connectors. (Sensor cable can be connected to sensors before or after placement). Different INVOS System sensors (adult, pediatric and neonatal) cannot be used on the same monitor.
5. Connect any optional accessories, including the Somanetics USB Flash Drive.




FOR DATA COLLECTION, THE USB FLASH DRIVE MUST BE CONNECTED TO THE INVOS SYSTEM BEFORE TURNING ON THE MONITOR. OR, THE USB FLASH DRIVE CAN BE CONNECTED ANYTIME, PROVIDED THE INVOS MONITOR IS RE-BOOTED TO INITIALIZE THE USB FLASH DRIVE BY PRESSING THE ON/OFF KEY. IF REMOVED DURING A CASE, THE ORIGINAL SOMANETICS USB FLASH DRIVE CAN BE RE-INSERTED TO CONTINUE DATA COLLECTION OF THE SAME CASE BY RE-BOOTING THE MONITOR AND PRESSING THE “PREVIOUS PATIENT” MENU OPTION ON THE START SCREEN.

6. Turn power ON by selecting the green  **ON/OFF** key.
7. After displaying the Welcome Screen, the INVOS System performs a 10-second self-test, stopping at the Start Screen.
8. If necessary, set the **DATE/TIME** by following the options menus. The date and time may only be changed from the Start Screen.
9. Press **NEW PATIENT**. Monitoring begins.
10. If desired, set alarms, scales, data channel labels and other modifiable preferences via the User Configuration Menu Screen.

Back Panel Connections and Elements



1. Digital Output Port **IOIOI**
2. VGA Output Port **IOI**
3. Potential Equalization \downarrow
(Grounding Post) Connector
4. Alarm Speaker and Cooling Fan
5. Fuse 

6. AC Mains (Master) Power Switch
7. AC Input, Power Cord Connector
8. Serial Number

Side Panel Connections and Elements




1. USB Connector 
2. Preamplifier A, Channel   Connector
3. Preamplifier B, Channel   Connector




To connect Preamplifier(s) to INVOS Monitor, align red dot on the silver connector with red dot on side panel input connection. Insert cable connector into side panel connection while keeping the dots aligned. Be sure to fully insert the connector until it locks.


To disconnect the Preamplifier, unlock the connection by drawing back on the ridged, silver outer sleeve of the connector. Then pull the cable out of the Preamplifier connection.



Start Screen: Date and Time Setting



- | | |
|--|---|
| 1. Software Version | 5. Date and Time on |
| 2. Navigation Bar (menu options) | Navigation Bar |
| 3. Blue  LED charge indicator | 6. Menu keys |
| 4. ON/OFF  key | 7. Home  key |

For changing the date or time, after the INVOS System is turned ON , follow these navigation menu options:

| | | | |
|-----------------|------------------|------------------|----------------------|
| NEW PATIENT | PREVIOUS PATIENT | DATE/TIME | NEXT MENU |
| DATE | TIME | | PREVIOUS MENU |
| DAY | MONTH | YEAR | PREVIOUS MENU |
| DAY = XX | INCREASE | DECREASE | PREVIOUS MENU |

Select DAY, MONTH or YEAR followed by the INCREASE or DECREASE options to adjust the date as necessary. Select HOURS, MINUTES or SECONDS followed by the INCREASE or DECREASE options to adjust the time as necessary. Press PREVIOUS MENU again to return to the Start Screen.



Cerebral Sensor Placement

1. To ensure good contact, clean/degrease the skin using mild soap and water. Ensure patient's skin is completely dry with a gauze pad.
2. Select sensor site on right and/or left sides of forehead (site selection will determine which region of the brain is monitored). Placement of the sensor in other locations, or over hair, may cause inaccurate readings, erratic readings, or no readings at all. Do not place the sensor over sinus cavities, the superior sagittal sinus, subdural or epidural hematomas or other anomalies such as arteriovenous malformations, as this may cause readings that are not reflective of brain tissue or no readings at all.
3. Peel off backing label from sensor and apply to the forehead as illustrated.
4. Secure the cable to a fixed object to avoid strain on the sensor to skin interface. (Sensor cable can be connected to sensors before or after placement. Hold blue connectors, not the black flex circuit.)
5. Press NEW PATIENT.

FOR EXTENDED MONITORING, SOMANETICS RECOMMENDS USING A NEW SOMASENSOR EVERY 24 HOURS AND FOLLOWING YOUR INSTITUTION'S GUIDELINES FOR SKIN INTEGRITY.

CAUTION: TO AVOID PRESSURE SORES, DO NOT APPLY EXTERNAL PRESSURE (E.G. HEADBANDS, WRAPS, TAPE) TO SENSOR.

SENSORS MAY OVERLAP, PROVIDED THE THREE OPTICAL WINDOWS ON EACH PAD REMAIN UNOBSTRUCTED.



Somatic Sensor Placement

Select sensor site over skeletal muscle tissue (site selection will determine which body region is monitored). Avoid placing the sensor over thick fatty deposits, hair or bony protuberances. Do not place the sensor over nevi, hematomas or broken skin, as this may cause readings that are not reflective of tissue or no readings at all.

Typical monitoring locations include:

- Renal area: over the latissimus dorsi muscle, T10-L2, posterior flank, right or left of midline
- Abdomen/gut: over the external abdominal oblique muscle
- Chest: over the pectoralis major muscle
- Forearm: over the brachioradialis muscle
- Upper arm: over the biceps muscle
- Calf: over the gastrocnemius muscle
- Upper leg: over the quadriceps muscle

Somatic Sensor Placement (continued)

FOR EXTENDED MONITORING, SOMANETICS RECOMMENDS USING A NEW SOMASENSOR EVERY 24 HOURS AND FOLLOWING YOUR INSTITUTION'S GUIDELINES FOR SKIN INTEGRITY.

CAUTION: TO AVOID PRESSURE SORES, DO NOT APPLY EXTERNAL PRESSURE (E.G. WRAPS, TAPE) TO SOMASENSOR.

SENSORS MAY OVERLAP, PROVIDED THE THREE OPTICAL WINDOWS ON EACH PAD REMAIN UNOBSTRUCTED.

1. To ensure good contact, wash area with mild soap and water and dry thoroughly.
2. Peel off backing and place SomaSensors. (Sensor cable can be connected to sensors before or after placement. Hold blue connectors, not the black flex circuit.)
3. Press NEW PATIENT.



User Configuration Menu Settings

The User Configuration Menu can be used to adjust a variety of settings including 2 or 4 channel monitoring, Alarms, rSO₂ Scale, Time Scale, USB Storage Rate, Event Mark Lists, color coding and naming of Channel Labels.

Assuming you have connected INVOS System components, the Monitor is ON and in Run Mode, press the HOME (🏠) key to display the Main Screen and follow these navigation bar options:

| | | | |
|---------------|---------------------------|--------------------|------------------|
| BASELINE MENU | EVENT MARK | ALARM AUDIO ON/OFF | NEXT MENU |
| OUTPUT SELECT | USER CONFIGURATION | TIME SCALE | PREVIOUS MENU |
| NEXT | PREVIOUS | INCREASE | DECREASE |

In the User Configuration Menu, select NEXT and PREVIOUS to scroll to the setting you wish to change. Adjust the setting.

When finished, press the HOME (🏠) key to return to the Main Screen. On the Main Screen, adjustments will appear.



User Configuration Menu Settings (continued)



Page Two

The User Configuration Menu default settings are:

| | DEFAULT | RANGE |
|------------------------------------|---------------------------------------|---|
| Upper Alarm Limit | OFF | 20-95 |
| Lower Alarm Limit | 40 | 15-90 |
| Alarm Volume | Medium | Low, Medium, High |
| Upper rSO ₂ Scale Limit | 100 | 60-100 in increments of 10; 5 above High Alarm Limit |
| Lower rSO ₂ Scale Limit | 30 | 0-30 in increments of 10; 5 below Low Alarm Limit |
| Trend (Time) Scale | 1 hour | 1, 2, 4, 8, 12, 24 hours |
| USB Storage Rate | 5 seconds | 5, 30, and 60 seconds |
| Event Mark List | OR | OR, ICU, VASCULAR, NICU |
| Channel Labels | L, R, S ₁ , S ₂ | A-Z, S ₁ -S ₄ |
| Color Coded Labels | NO | YES, NO |

User Configuration Menu (2- or 4-Channel Monitoring)

THE INVOS SYSTEM MUST BE RESTARTED FOR THIS CHANGE TO TAKE EFFECT.

To change the software version from two to four channel monitoring (assuming you have connected INVOS System components, the Monitor is ON and in Run Mode), press the HOME (🏠) key to display the Main Screen and follow these navigation bar options:

| | | | |
|---------------|---------------------------|--------------------|------------------|
| BASELINE MENU | EVENT MARK | ALARM AUDIO ON/OFF | NEXT MENU |
| OUTPUT SELECT | USER CONFIGURATION | TIME SCALE | PREVIOUS MENU |
| NEXT | PREVIOUS | INCREASE | DECREASE |
| NEXT | PREVIOUS | INCREASE | DECREASE |
| NEXT | PREVIOUS | YES | NO |
| NEXT | PREVIOUS | INCREASE | DECREASE |

In the User Configuration Menu, select NEXT to scroll to the 2 OR 4 CHANNELS setting. Select INCREASE and DECREASE to choose the number of channels monitored and displayed.

When finished, press the HOME (🏠) key to return to the Main Screen. **Press the ON/OFF (🔌) key twice to restart the system.** When the Monitor has returned to Run Mode, the number of channels selected should be displayed.

SEPARATE SETS OF USER CONFIGURATION PARAMETERS ARE STORED FOR 2- AND 4-CHANNEL MONITORING. PREVIOUS SETTINGS FOR 2-CHANNEL MONITORING WILL RESUME WHEN THIS CONFIGURATION IS SELECTED, AND VICE VERSA.

Main Screen (Two-Channel Monitoring)



1. Date and Time
2. Adult, Pediatric or Neonatal Indication
3. Upper and Lower Alarm Limits
4. Trend Data Graph (Channel 1)
5. Trend Data Graph (Channel 2)
6. rSO₂ Trend Scale (High & Low)
7. Time Scale; in hours
8. Menu Options (Navigation Bar)
9. Alarm Status Message
10. Alarm Symbol
11. Signal Strength Indicator
12. rSO₂ Baseline and Symbol
13. rSO₂ Relative (%) change from Baseline
14. Real-Time rSO₂ (index) Values
15. Data Channel Label (e.g. L; left cerebral, R; right cerebral, S; somatic)
16. Event Mark line and numerical code
17. Battery Power Status

Main Screen (Four-Channel Monitoring)



- | | |
|--|---|
| 1. Date and Time | 10. Alarm Symbol |
| 2. Adult, Pediatric or Neonatal Indication | 11. Signal Strength Indicator |
| 3. Upper and Lower Alarm Limits | 12. rSO ₂ Baseline and Symbol |
| 4. Trend Data Graph (for Channels 1 & 2) | 13. rSO ₂ Relative (%) change from Baseline |
| 5. Trend Data Graph (for Channels 3 & 4) | 14. Real-Time rSO ₂ (index) Values |
| 6. rSO ₂ Trend Scale (High & Low) | 15. Data Channel Label (e.g. L; left cerebral, R; right cerebral, S; somatic) |
| 7. Time Scale; in hours | 16. Event Mark line and numerical code |
| 8. Menu Options (Navigation Bar) | 17. Battery Power Status |
| 9. Alarm Status Message | |

B Setting rSO₂ Baselines

An rSO₂ Baseline should be set so subsequent rSO₂ changes from baseline can be measured. rSO₂ numbers will turn red if rSO₂ is above or below set alarm thresholds. Baseline status will turn red if rSO₂ is ≤ 20% below set baseline.

In the Baseline Menu, baseline(s) may be set for all channels in use at once, individually for selected channel(s) or manually to a previous rSO₂ value.

IF A BASELINE CANNOT BE SET PRE-INDUCTION FOR PEDIATRIC PATIENTS, IT IS RECOMMENDED THAT A VALUE OF 50 IS USED FOR THE DESATURATION THRESHOLD.

Setting an Awake rSO₂ Baseline for All Channels: Recommend Prior To Pre-op or O₂ Mask

As soon as rSO₂ trend data begins to be gathered and patient condition is stable, an awake baseline reading should be obtained (e.g. for surgical patients, prior to induction).

To set a baseline for all channels, press the HOME (🏠) key to display the Main Screen, then follow these navigation menu options:

| | | | |
|----------------------|-------------------|--------------------|------------|
| BASELINE MENU | EVENT MARK | ALARM AUDIO ON/OFF | NEXT MENU |
| SET BASELINES | RESTORE BASELINES | SET CHANNEL | MANUAL SET |

Set Individual Channels to Current rSO₂ Baseline

To set an individual channel baseline, press the HOME (🏠) key to display the Main Screen, then follow these navigation menu options:


| | | | |
|----------------------|----------------------|----------------------|----------------------|
| BASELINE MENU | EVENT MARK | ALARM AUDIO ON/OFF | NEXT MENU |
| SET BASELINES | RESTORE BASELINES | SET CHANNEL | MANUAL SET |
| CHANNEL 1 SET | CHANNEL 2 SET | CHANNEL 3 SET | CHANNEL 4 SET |

Select the channel(s) to set the baseline at the current rSO₂ value. This will display an event mark line and code on the display screen. When finished, press the HOME (🏠) key to exit the Baseline Menu and return to Main Screen.

Setting rSO₂ Baselines (continued)

WHEN USING THE MANUAL SET BASELINE FEATURE, THE NEW BASELINE VALUES WILL NOT BE REFLECTED IN ANY REAL-TIME DATA OUTPUT, INCLUDING TO THE USB FLASH DRIVE, OR WHEN REVIEWING CASES IN CASE ARCHIVE MODE. FOR ARCHIVING PURPOSES, THE CASE SHOULD BE DOWNLOADED IN CASE HISTORY, BEFORE TURNING OFF THE MONITOR, IMMEDIATELY AT THE COMPLETION OF THE CASE.

Manually Set Baselines to Previous rSO₂ Value

The manually set baseline feature is used when previous rSO₂ values are to be set as baselines (e.g., an awake baseline was forgotten to be set prior to induction). To select and manually set previous rSO₂ values, press the HOME  key to display the Main Screen, then follow these navigation menu options:

| | | | |
|----------------------|----------------------|--------------------|-------------------|
| BASELINE MENU | EVENT MARK | ALARM AUDIO ON/OFF | NEXT MENU |
| SET BASELINES | RESTORE BASELINES | SET CHANNEL | MANUAL SET |
| NEXT PAGE | PREVIOUS PAGE | LAST PAGE | SCROLL MENU |

Select **NEXT PAGE**, **PREVIOUS PAGE** or **LAST PAGE** to review pages of rSO₂ data and locate the desired time and value to set baselines to.


Once the time and data on the page being displayed has been located, follow these navigation bar menu options:


| | | | |
|-------------|-----------------|-------------------------|--------------------|
| NEXT PAGE | PREVIOUS PAGE | LAST PAGE | SCROLL MENU |
| NEXT | PREVIOUS | SELECT BASELINES | PREVIOUS MENU |

Select **NEXT** or **PREVIOUS** to scroll to the desired line of rSO₂ data, then **SELECT BASELINES** to manually set the baseline.

Setting rSO₂ Baselines (continued)

Restoring a Baseline/Replacing a Sensor

To restore a baseline or replace a sensor, press the HOME  key to display the Main Screen. If necessary, select the ALARM AUDIO ON/OFF option to deactivate audible alarms.

CAUTION: SELECTING THE ALARM AUDIO ON/OFF  OPTION WILL SILENCE THE AUDIBLE ALARMS PERMANENTLY UNTIL RE-ACTIVATED OR THE INVOS SYSTEM IS REBOOTED.

FOR CONTINUOUS MONITORING WHILE CHANGING MULTIPLE SENSORS, REMOVE ONE SENSOR FOLLOWED BY SENSOR PREP AND PLACEMENT INSTRUCTIONS BEFORE REPLACING OTHER SENSORS.


Prep the patient, apply sensors, then follow these navigation menu options:

| | | | |
|----------------------|--------------------------|--------------------|------------|
| BASELINE MENU | EVENT MARK | ALARM AUDIO ON/OFF | NEXT MENU |
| SET BASELINES | RESTORE BASELINES | SET CHANNEL | MANUAL SET |


Select ALARM AUDIO ON/OFF to re-activate audible alarms, if necessary.

Event Marking

An Event Mark may be used to mark significant occurrences. It is displayed on the screen or stored in memory for later review. A vertical dashed line and event code will appear on the trend graph at the marked time. An event code will also appear at that time in the line of tabular trends screen or stored and outputted data.

To mark an event, press the HOME  key to display the Main Screen, then follow these navigation menu options:

| | | | |
|-------------------|-----------------------|---------------------|---------------|
| BASELINE MENU | EVENT MARK | ALARM AUDIO ON/OFF | NEXT MENU |
| NEXT EVENT | PREVIOUS EVENT | SELECT EVENT | PREVIOUS MENU |

From the Event Mark List, select NEXT EVENT and PREVIOUS EVENT to highlight the desired event code or description. SELECT EVENT to choose the highlighted event code and return to the trend data graph. Press the HOME  key to return to the trend data graph without storing an event code.

OR EVENT MARK LIST

- | | | | |
|----|-----------------------|----|---------------------------|
| 1 | Miscellaneous | 21 | Afterload Reduction |
| 2 | Set baseline | 22 | Blood Transfusion |
| 3 | Intubated | 23 | Cardioversion |
| 4 | Sternotomy | 24 | Cell saver blood |
| 5 | Cannulate | 25 | Cerebral Perfusion On |
| 6 | On CPB | 26 | ECLS On |
| 7 | Cross Clamp On | 27 | FFP / Platelets |
| 8 | Cooling | 28 | Fluid/Volume Expander |
| 9 | Cardioplegia | 29 | Hemoconcentrate / MUF |
| 10 | Warming | 30 | Inotrope |
| 11 | Cross Clamp Off | 31 | Increase Anesthetic |
| 12 | Off CPB | 32 | Increase CO ₂ |
| 13 | Closing Sternum | 33 | Increase FiO ₂ |
| 14 | Arrhythmia | 34 | Increase Pump Flow |
| 15 | Circulatory Arrest | 35 | Paced |
| 16 | Hypocapnia | 36 | Reposition Cannula |
| 17 | Hypotension | 37 | Reposition Clamp |
| 18 | One Lung Ventilation | 38 | Reposition Head |
| 19 | Pump Flow Down | 39 | Reposition Heart |
| 20 | Reduced Venous Return | 40 | Vasopressor |

Event List Key Code:

White events - reflect relatively routine occurrences during surgery. These are listed in the order they typically occur during a case.

Yellow events - mark events clinicians may consider “cautionary” because they could cause ischemia. These are listed alphabetically.

Green events - mark interventions taken and are listed alphabetically.

PEDIATRIC/ADULT ICU EVENT MARK LIST

| | |
|------------------------|------------------------|
| 41 Miscellaneous | 61 Blood Transfusion |
| 42 Set Baseline | 62 Chest Closed |
| 43 Enteral feeding | 63 Dialysis / CRRT |
| 44 Extubated | 64 Diuretic |
| 45 Intubated | 65 ECLS On |
| 46 Reposition Patient | 66 ECLS Circuit Change |
| 47 Sensor Change | 67 ECLS Off |
| 48 Apnea | 68 ET Tube Suctioned |
| 49 Arrhythmia | 69 Fluid Bolus |
| 50 Bradycardia | 70 FFP/Platelets |
| 51 Cardiac Arrest | 71 Hi Frequency Vent |
| 52 ICP Changes | 72 Hypothermia |
| 53 LOC Changes | 73 Inotrope |
| 54 Painful Procedure | 74 Nitric Oxide |
| 55 Seizure Activity | 75 Paralytic |
| 56 Tamponade | 76 PDA Ligated |
| 57 Afterload Reduction | 77 Prostaglandin |
| 58 Anti-Arrhythmic | 78 Sedation |
| 59 Anti-Epileptic | 79 Vasopressor |
| 60 Anti-Pyretic | 80 Ventilator Change |

Event List Key Code:

White events - reflect relatively routine occurrences. These are listed in the order they typically occur.

Yellow events - mark events clinicians may consider "cautionary" because they could cause ischemia. These are listed alphabetically.

Green events - mark interventions taken and are listed alphabetically.

NICU EVENT MARK LIST

| | | | |
|-----|-----------------------|-----|-----------------------|
| 121 | Miscellaneous | 141 | Cooling Cap On-Off |
| 122 | Set Baseline | 142 | Dialysis/CRRT |
| 123 | Enteral Feeding | 143 | ECLS On |
| 124 | Extubated-Intubated | 144 | ECLS Circuit Change |
| 125 | Bag Mask Ventilation | 145 | ECLS Off |
| 126 | Conventional Vent | 146 | Fluid Bolus |
| 127 | Hi Frequency Vent | 147 | Fem Art CutDwn Start |
| 128 | Position Change | 148 | Fem Art CutDwn Stop |
| 129 | Sensor Change | 149 | Hypothermia |
| 130 | Weigh Patient | 150 | Nitric Oxide On – Off |
| 131 | Apnea | 151 | NG Tube In – Out |
| 132 | Arrhythmia | 152 | Paralytic |
| 133 | Cardiac Arrest/CPR | 153 | PDA Ligated |
| 134 | ICP Changes | 154 | Prostaglandin |
| 135 | LOC Changes | 155 | OR Procedure Bedside |
| 136 | Seizure | 156 | Sedation |
| 137 | Tamponade | 157 | Suction ET Tube |
| 138 | Anti-Arrhythmic | 158 | Vasopressor |
| 139 | Anti-Epileptic | 159 | Ventilator Change/BMV |
| 140 | Blood Trans/Platelets | 160 | Whole Body Cooling |

Event List Key Code:

White events - reflect relatively routine occurrences. These are listed in the order they typically occur.

Yellow events - mark events clinicians may consider “cautionary” because they could cause ischemia. These are listed alphabetically.

Green events - mark interventions taken and are listed alphabetically.

Alarms and Status Messages with rSO_2 Values Displayed

During normal monitoring situations when rSO_2 values are being displayed, alarm and status messages may appear at the top of the Main Screen. Additional operating indications to assist the user will appear on the display screen at the right of the affected channel label (e.g. 2-digit rSO_2 values turn red when value is below baseline desaturation threshold.)

| STATUS MESSAGE | COLOR, LOCATION | CAUSE | ACTION |
|---|--|--|---|
| ALARM HIGH <i>accompanied by alarm tone</i> STATUS CODE: 5 | Red, top right | Upper Alarm Limit is exceeded. | Profile Patient |
| ALARM LOW <i>accompanied by alarm tone</i> STATUS CODE: 6 | Red, top right | Lower Alarm Limit is exceeded. | Profile Patient |
| BATTERY LOW <i>accompanied by alarm tone</i> STATUS CODE: 9 | Red, bottom right above navigation bar | Battery power is critical. | Plug in A/C Power and turn master (mains) switch on. |
| COM PORT UNAVAILABLE STATUS CODE: 12 | Yellow, top right | Output device not connected. | Check connections with Com Port device (PC). |
| USB FULL STATUS CODE: 20 | Yellow, top right | USB Flash Drive memory is full. | Insert a blank Flash Drive. |
| USB STORAGE ERROR STATUS CODE: 14 | Yellow, top right | USB Flash Drive error. | Check USB Flash Drive connection. Insert a Flash Drive. |
| DATE AND TIME | Green, top left | | |
| ADULT, PEDIATRIC or NEONATAL | Green, top left | Adult, Pediatric or Neonatal application dependent on sensor model used. | |
| SYSTEM SIGNAL OK STATUS CODE: 4 | No text message displayed on screen | SSI at full strength | |

Alarms and Status Messages without rSO₂ Values Displayed

During situations when rSO₂ values cannot be calculated, status messages may appear in place of the rSO₂ values at the right of the corresponding affected channel label indicator.

| STATUS MESSAGE | CAUSE | ACTION |
|--|--|---|
| SENSOR NOT CONNECTED <i>accompanied by alarm tone</i> STATUS CODE: 1 | Disconnect at sensor (patient) or cable. No rSO ₂ value will appear. | Check Sensor to Reusable Sensor Cable connection. Check Sensor Cable connection to the Preamplifier. Restore baseline values. |
| EXCESSIVE LIGHT <i>accompanied by alarm tone</i> STATUS CODE: 2 | Sensor has lifted off the patient, or there is too much outside light. No rSO ₂ will appear. | Check the connection to patient and replace sensor, if necessary. If the sensor is properly adhered to the patient, try loosely draping the sensor. |
| POOR SIGNAL QUALITY <i>accompanied by alarm tone</i> STATUS CODE: 3 | INVOS System is operating, but rSO ₂ values will not appear because they are unstable or may be corrupted by a noisy power source or a very weak signal. | Try another hospital grade electrical outlet. Relocate sensor, if possible. |
| PREAMP NOT CONNECTED STATUS CODE: 11 | Preamplifier is not connected to the INVOS Monitor. | Make sure the Preamplifier connector is locked into the side of the monitor. Normal operation of the INVOS System will resume when the Preamplifier is reconnected. |
| REPLACE SENSOR <i>accompanied by alarm tone</i> STATUS CODE: 17 | INVOS System cannot read calibration data from the sensor. Mixed adult, pediatric and neonatal sensor models, moisture intrusion, a defective sensor, or a defective Reusable Sensor Cable could cause this. | Check for correct sensor model type. Connect new SomaSensor. Replace Reusable Sensor Cable. |
| INTERFERENCE DETECTED <i>accompanied by alarm tone</i> STATUS CODE: 19 | Noise corrupting rSO ₂ data (e.g. electrocautery). | Move or turn off source of interference. Normal operation of the INVOS System will resume when the excessive noise stops. |

Troubleshooting Signal Strength Indicator (SSI)

The Signal Strength Indicator (SSI) reflects the strength and stability of the signal that is used to make an rSO₂ measurement. The SSI depends on physiological factors (different from individual to individual) and possible disturbances from secondary sources/influences in the monitor's environment. An inherent benefit of the technology is the ability of the device to determine rSO₂ values from very small signals.


The chart below provides SSI ranges and a list of the possible causes and action to correct for each.

| SSI INDICATION | DESCRIPTION/CAUSE | ACTION |
|---|---|---|
| 5 bars <i>STATUS CODE: 4</i> | Full strength signal. | No action required. |
| 1 through 4 bars - stable | Signal is weaker, but is strong and stable enough to generate an accurate rSO ₂ value. | No action required. |
| 1through 4 bars - fluctuating and unstable | Signal is weak and unstable. | |
| | Signal Interference. | Identify any possible sources of interference, such as electro-cautery or similar electrical equipment interferences, intense light or radiation sources, moving artifacts or perturbations on the mains. |
| | Sensor application over nevi, sinus cavities, the superior sagittal sinus, subdural or epidural hematomas or other anomalies such as arteriovenous malformations, broken skin, thick fatty deposits, hair or bony protuberances, or area contaminated with residue. | Clean the area, re-apply SomaSensor to patient in an adjacent area; follow the sensor application instructions and make sure no hair, anomalies or residues are under the sensor. |
| | Defective Sensor. | Apply new Sensor to patient. |
| | Moisture in sensor connection. | Dry the cable and sensor connectors. Make sure the connection is not exposed to moisture. |
| 0 bar | See Error Messages in data channel area on display screen (See Status Messages – without rSO ₂ Values Displayed) | Follow on screen instructions. |

Reviewing Previous Cases on the Monitor

The INVOS System offers several different options to store, output and review data collected. One option is Case Archive Mode, allowing review of up to 28 cases of 24-hour duration. When data storage exceeds 24-hours or midnight is reached, a new case will be started. Case History data files are selected for review in the Case Archive File List.


Data files created by the INVOS System are automatically named by the software using the date (YYMMDD) followed by a letter indicating the type of recording (N = Normal Run Mode, C = Control, I = Intervention) and a file extension consisting of a letter (H = History, or R = Real-time) and a sequence number. For example, a CASE HISTORY file recorded on November 13, 2007 in Normal Run Mode would have a filename of 071113N.H1. An asterisk (*) indicates a desaturation of 20% or more from baseline or ≤ 40 rSO₂ values occurred during the case.



Reviewing Previous Cases on the Monitor


(continued)

MORE THAN ONE (1) HOUR OF DATA MUST BE STORED TO SCROLL IN CASE ARCHIVE MODE.

To access Case Archive Mode in graphical format, press the ON/OFF  key to display the Start Screen and follow these navigation menu options:

| | | | |
|---------------------|----------------------|--------------------|------------------|
| NEW PATIENT | PREVIOUS PATIENT | DATE/TIME | NEXT MENU |
| CASE ARCHIVE | PREVIOUS MENU | | |
| NEXT FILE | PREVIOUS FILE | SELECT FILE | START SCREEN |
| ←TREND GRAPH | TREND GRAPH→ | FULL GRAPH | NEXT MENU |

Select NEXT FILE or PREVIOUS FILE to scroll through the Case Archive File List. SELECT FILE will load and display the entire patient file graphically. Select ← TREND GRAPH and TREND GRAPH → to scroll in one-hour increments. Select FULL GRAPH to display the entire case.




Downloading Cases Via USB Flash Drive

The USB Flash Drive allows data from the INVOS System to be collected and transferred to other devices such as a PC. Data may be collected in REAL-TIME or CASE HISTORY formats.


Once the INVOS System and USB Flash Drive are connected and in Run Mode, the Monitor will begin data outputting to the USB Flash Drive automatically.



Up to 24 hours of data may be stored per case file. When data storage exceeds 24 hours or midnight is reached, a new case file will be started.

FOR DATA COLLECTION, THE USB FLASH DRIVE MUST BE CONNECTED TO THE INVOS SYSTEM BEFORE TURNING ON THE MONITOR. OR, THE USB FLASH DRIVE CAN BE CONNECTED ANYTIME, PROVIDED THE INVOS MONITOR IS RE-BOOTED TO INITIALIZE THE USB FLASH DRIVE BY PRESSING THE ON/OFF  KEY.

The stored data can be accessed and graphed on a personal computer using a common spreadsheet program like Microsoft Excel. The file will be stored as an ASCII text file. See the Operations Manual for the data format.


Downloading an Individual Case from Archive List

Access Case Archive Mode to download a desired file in the Case Archive File list. Assuming you have connected INVOS System components and connected the USB Flash Drive, press the ON/OFF  key to display the Start Screen and follow these navigation menu options:

| | | | |
|--|---|--------------------|------------------|
| NEW PATIENT | PREVIOUS PATIENT | DATE/TIME | NEXT MENU |
| CASE ARCHIVE | START SCREEN | | PREVIOUS MENU |
| NEXT FILE | PREVIOUS FILE | SELECT FILE | NEXT MENU |
|  TREND GRAPH | TREND GRAPH  | FULL GRAPH | NEXT MENU |
| TABULAR TRENDS | DOWNLOAD | PREVIOUS MENU | NEXT MENU |
| USB STORAGE | DIGITAL OUTPUT | | PREVIOUS MENU |

Downloading Cases Via USB Flash Drive (continued)


Downloading All Stored Cases from Archive List

Access Case Archive Mode to download all files in the Case Archive File list. Assuming you have connected INVOS System components and connected the USB Flash Drive, press the ON/OFF  key to display the Start Screen and follow these navigation menu options:

| | | | |
|---------------------------|-------------------|---------------|------------------|
| NEW PATIENT | PREVIOUS PATIENT | DATE/TIME | NEXT MENU |
| CASE ARCHIVE | START SCREEN | | PREVIOUS MENU |
| NEXT FILE | PREVIOUS FILE | SELECT FILE | NEXT MENU |
| DOWNLOAD FILES | DELETE FILES | START SCREEN | PREVIOUS MENU |
| DOWNLOAD ALL FILES | CASE ARCHIVE MENU | PREVIOUS MENU | |

All files except EXAMPLE.1 and EXAMPLE.2 will be transferred to the USB Flash Drive in a maximum of four minutes (28 cases of 24-hour duration). As the files are downloading, the file currently being transferred appears above the Case Archive File List.

Downloading Case Data at the End of a Case

Assuming you have connected INVOS System components and connected the USB Flash Drive, press the HOME  key to display the Start Screen and follow these navigation menu options:

| | | | |
|----------------------|--------------------|--------------------|------------------|
| BASELINE MENU | EVENT MARK | ALARM AUDIO ON/OFF | NEXT MENU |
| OUTPUT SELECT | USER CONFIGURATION | | PREVIOUS MENU |
| USB | DIGITAL OUTPUT | REVIEW | PREVIOUS MENU |
| CASE HISTORY | | PREVIOUS MENU | MAIN MENU |


The Case History data will output until complete.

SOMANETICS®

Site-Specific Saturation Data With Near-infrared Spectroscopy (NIRS)



Four data channels enable simultaneous monitoring of cerebral and somatic regional oxygen saturation (rSO_2)



The INVOS System uses two depths of light penetration to subtract out surface data, resulting in a regional oxygenation value for deeper tissues.

Regional Oximetry Vs. Other Oximetry

Regional (Capillary) Oximetry (rSO₂) Clinical Characteristics

- Non-invasive
- Capillary (venous and arterial) sample
- Measures the balance between O₂ supply and demand
- End organ oxygenation and perfusion
- Requires neither pulsatility nor flow

Pulse (Arterial) Oximetry (SpO₂) Clinical Characteristics

- Noninvasive
- Arterial sample
- Measures O₂ supply in the periphery
- Systemic oxygenation
- Requires pulsatility and flow

Central (Venous) Oximetry (SvO₂) Clinical Characteristics

- Invasive
- Venous sample
- Measures O₂ surplus in central circulation
- Systemic oxygen reserve
- Requires flow


Factors That Can Affect Cerebral Oxygenation

Tailoring Oxygen Delivery

- Mean Arterial Pressure - rSO_2 values help to tailor individual treatment.
- CO_2 - Higher CO_2 causes cerebral vasodilation to allow greater blood flow to the brain.
- Cardiac Index and Pump Flow - Insufficient CI and pump flow cause rSO_2 declines.
- Blood Products - Hemodilution has a negative effect on oxygen carrying capacity; the addition of red blood cells increases rSO_2 .
- Mechanical Issues - Cannula misplacement, neck position and other issues can cause rapid and dramatic rSO_2 changes.

Tailoring Oxygen Consumption

- Anesthetic Agent and Depth - Anesthesia decreases O_2 consumption by the brain; light anesthesia may require more O_2 delivery.
- Temperature - Temperature changes can cause rapid rSO_2 changes. Cooling usually increases rSO_2 , while warming usually lowers it.



rSO₂ Reflects Oxygen Balance

rSO₂

- Increases with rise in delivery or fall in demand
- Decreases when delivery falls or if there is an uncompensated rise in demand

Oxygen Delivery/Supply

Influenced by:

- Oxygen Content
 - Hemoglobin concentration
 - Hemoglobin saturation
- Cardiac Output
 - Optimize heart rate
 - Idealize preload
 - Improve contractility
 - Manipulate afterload

Oxygen Demand/Consumption

Increased by:

- Fever, shivering
- Malignancy, severe infection
- Cold stress (neonates)
- Seizures, status epilepticus
- Wounds and burns
- Pain

Decreased by:

- Hypothermia, without shivering
- Sedation and paralysis
- Shunting or decreased extraction

The Cerebral - Somatic Relationship

The INVOS System provides perfusion data from vascular beds that represent opposite poles of regional circulation and have different extraction ratios.



- Cerebral
 - High flow, high extraction organ
 - Compensatory mechanisms
 - autoregulation
 - flow-metabolism coupling
 - Cerebral desaturations are a *late* indicator of shock if cerebral autoregulation is intact



- Somatic
 - Variable flow, lower O₂ extraction
 - Flow is highly influenced by autonomic (sympathetic) tone
 - Somatic desaturations may be an *early* indicator of shock (i.e., peripheral circulation is shutting down to preserve the brain)

rSO₂ Targets and Thresholds

Targets and thresholds are expressed in rSO₂ numerical values and % changes from baseline. This enables customized management based on each patient's unique physiology and clinical anomalies. As such, each patient serves as his/her own control rather than being relegated to a fixed protocol.


Cerebral

- High flow, high extraction organ
- Typical range: 60-80
- ≤ 50 rSO₂ or 20% Δ from rSO₂ baseline is common intervention trigger
- ≤ 45 rSO₂ or 25% Δ from rSO₂ baseline is critical threshold

Somatic

- Variable flow, lower O₂ extraction
- 5-20 points higher than cerebral rSO₂; changes in this cerebral - somatic variance may be indicative of pathology.

rSO₂ Reactions to Common Intra-Op Events



This aortic arch reconstruction with deep hypothermic circulatory arrest demonstrates common times when cerebral perfusion can be impacted both positively and negatively. Adequate rSO₂ levels can provide confidence that the brain is well perfused, while dips in rSO₂ near or past the patient's desaturation threshold indicate opportunities for therapeutic interventions.

Operating Room Interventions to Improve rSO_2

While each hospital will have its own care protocols, the following interventions have been shown to improve rSO_2 .

Perfusion imbalance

- Blood pressure
- Mechanical obstruction (cannula position)
- Increase cardiac output (pump flow)
- Plasma volume expansion
- CO_2 to physiologic level

Limited ischemia tolerance


- Neuroprotective Agent
- Additional cooling

Dysoxygenation

- Increase FiO_2
- Whole blood
- Reintubate

Insufficient Anesthetic Depth

- Increase anesthetic delivery



PICU Interventions to Improve rSO_2 - Cerebral

While each hospital will have its own care protocols, the following interventions have been shown to improve rSO_2 .

Cerebral

Increase cerebral perfusion pressure

- Increase arterial blood pressure
- Increase systemic vascular resistance
- Increase cardiac output
- Reduce central venous pressure

Increase arterial oxygen content


- Transfuse red cells
- Raise arterial partial pressure of oxygen

Reduce cerebral vascular resistance

- Raise arterial partial pressure of carbon dioxide

Reduce cerebral metabolic rate

- Control hyperthermia
- Sedation



PICU Interventions to Improve rSO_2 - Somatic

While each hospital will have its own care protocols, the following interventions have been shown to improve rSO_2 .

Somatic

Increase total cardiac output

- optimize: preload, afterload, rate, rhythm and contractility


Reduce sympathetic outflow

- increase inotropes
- decrease Nitroprusside

Increase hematocrit

Maintain temperature WNL

Consider regional vasodilation



Patient Care

The care team should follow its institution's own protocols for patient monitoring and skin integrity. Some examples pertaining to INVOS System use include:

- Change sensors routinely
 - Recommended every 24 hours
 - Assess skin integrity
 - Skin under and around the sensor
 - Consider allowing the skin to “breathe”
 - Change one sensor at a time to ensure continuous monitoring
 - Place new sensor in the same location for continuity of measurements

Charting

The care team should follow its institution's own protocols for patient charting. Some examples pertaining to INVOS System use include:

- Report and record rSO_2 desaturations below 50, or below a 20% drop from the patient's baseline
- Report and record interventions taken, and your evaluation of interventions
- Share current and past rSO_2 trend data at shift changes, physician rounds, etc.
- Collaborate on appropriate treatment options




Interfering Factors

rSO₂ values can be affected by:

- Jaundice (hyper-bilirubinemia)
- Severe tissue edema
- Major sinus problems (cerebral)

Reversal of Shock




This infant with septic/cardiogenic shock and severe global hypoperfusion presented with below critical rSO₂ levels. Subsequent rSO₂ values reflect patient responses to a series of interventions. Regional oxygenation was normalized at higher levels, including restoration of typical somatic oxygenation at 5 - 20 points greater than cerebral.¹⁻³

¹ Hoffman GM, et al. Semin Thorac Cardiovasc Surg Pediatr Card Surg Annu. 2005;12-21.

² Petrova A and Mehta R. Pediatr Crit Care Med 2006;7:449-454.

³Data on file.


Detection of Mechanical Obstructions - 5 kg Infant



During this BT shunt case, initial intubation was difficult and the ET tube later became dislodged. This mechanical obstruction caused an immediate ~20 point drop in cerebral rSO₂ to a critical level of 45. SpO₂ values remained normal until about 1 minute later.


Similarly, after cross clamping the right subclavian artery, cerebral rSO₂ fell ~40 points over a period of 12 minutes, again dipping below critical. The perfusionist changed flow and pressure with no improvement. The surgeon then checked the cross clamp which had inadvertently included the innominate artery occluding some flow to the brain. Upon clamp repositioning, cerebral rSO₂ immediately improved. There were no other hemodynamic indicators that the innominate artery was partially occluded. Without cerebral rSO₂ values, the brain would have had inadequate oxygenation for ~25 minutes.

Effect of Organ Compression / Manipulation - Two-day-old Infant



This giant omphalocele repair case emphasizes the potentially negative effects of organ compression and manipulation on regional oxygen saturation. Somatic monitoring enabled clinicians to visualize compromise to peri-renal perfusion in real-time and adjust their surgical approach to provide immediate relief of abdominal compression. Interventions included repositioning the omphalocele, spacing their manipulations to allow for perfusion recovery in between, vent changes and PRBCs to replace losses. Without placement of the peri-renal sensor, prolonged desaturation of the kidney area may have continued without notice for a significant period of time, increasing the risk of post-op renal dysfunction. This same desaturation threat can occur in similar compression-prone scenarios such as lifting/manipulating the heart and closure of the chest.

Oscillations and Desaturation Associated with Seizures



1-5. Episodes of Seizure; data recorded every 10 seconds for best resolution

An EEG monitoring study detected seizure activity in 26% of circulatory arrest infants during and after surgery, which can result in serious developmental consequences.¹ This cerebral oximetry trend graph is marked at the times of EEG documented seizures.² Note the transition from a smooth line to one with oscillations and a desaturation. The INVOS System can help identify potential seizure activity--even silent seizures and hemisphere of origin.³⁻⁴

¹ Newburger JW, et al. *New Engl J Med* 1993;329:1057-1064.

² Courtesy of HL Edmonds, Jr, PhD


³ Diaz GA, et al. *Eur J Paediatr Neurol*. 2006;10(1):19-21. Epub 2006 Mar 10.

⁴ Shuhaiber H, et al. *J Child Neurol*. 2004 Jul;19(7):539-40.


ECMO Management

rSO₂ values are useful in determining the need to initiate ECMO therapy as well as in providing immediate patient response to separation from ECMO.


Initiating ECMO



Separation from ECMO



Sole Early Indicator of Low Cardiac Output Syndrome - 10 kg Child



Roughly 24 hours after ASD/ VSD repair, somatic and cerebral rSO₂ values frequently dipped to sub-optimal levels, while venous and pulse oximetry remained stable. Nursing and respiratory staff documented bilateral breath sounds throughout the night. By early morning, a critical desaturation of both SpO₂ and rSO₂ occurred with recovery of SpO₂ to only 88% despite increasing FiO₂ to 100%. A stat CXR was obtained revealing significant pleural effusion that was tapped and drained. This enabled the heart to pump sufficiently to perfuse the somatic/peri-renal area, as evidenced by the dramatic rise in somatic rSO₂. This case demonstrates the unique value of site-specific oxygenation data to aid in the detection of low cardiac output syndromes and limited regional perfusion.

SOMANETICS®

Ordering Information

| <i>Description</i> | <i>Model Number</i> |
|---|---------------------|
| Disposable Adult SomaSensor | SAFB-SM |
| Disposable Pediatric SomaSensor | SPFB |
| INVOS 5100C Cerebral/Somatic Oximeter | 5100C-USA |
| Preamplifier with Cable, Channel 1 & 2 (Preamp) | 5100C-PA |
| Preamplifier with Cable, Channel 3 & 4 (Preamp) | 5100C-PB |
| Reusable Sensor Cable, Channel 1 (Blue) | RS-1 |
| Reusable Sensor Cable, Channel 2 (Taupe) | RS-2 |
| Reusable Sensor Cable, Channel 3 (Orange) | RS-3 |
| Reusable Sensor Cable, Channel 4 (Green) | RS-4 |
| Power Cord | 5100-PCNAM |
| USB Flash Drive | 5100C-USB |
| VueLink Adapter Cable (for connecting to a Philips VueLink Interface Module) | VLI |
| Null Modem Cable (for digital output) | 312170 |
| Portable Roll Stand (for 5100C) | 5100C-RS |
| Swivel Arm (for 5100C) | 5100C-SA |
| Operations Manual | 5100C-M-USA |
| Quick Reference Guide for Adults | 5100C-RFA |
| Quick Reference Guide for Pediatrics | 5100C-RFP |

Accessories can be ordered by contacting Somanetics' Customer Service at 800.359.7662 in the U.S. or 248.689.3050, ext. 255 outside the United States or via the Internet at www.somanetics.com.

Warning: Accessories not supplied by Somanetics may not meet EN60601-1-2 (IEC 601-1-2) standards. Contact Somanetics' Customer Service Department for compatible products that may meet these requirements.



© Somanetics Corporation, Troy, Michigan. All rights reserved.

Somanetics, INVOS, SomaSensor and "Reflecting the color of life" are registered trademarks of Somanetics Corporation. NIRSensor is a trademark of Somanetics Corporation. All other trademarks that appear in these instructions are property of their respective owners.

U.S. patents 5139025, 5140989, 5217013, 5349961, 5465714, 5477853, 5482034, 5584296, 5697367, 5795292, 5902235, 6615065. Foreign patents. Other patents pending.

Rx ONLY

PN 313380, rev. A

Somanetics Corporation • 1653 East Maple Road, Troy, MI, U.S.A. 48083
Phone: 800.359.7662 • 248.689.3050, ext. 255 • www.somanetics.com